

TransMAX: Maximizing Adenovirus uptake



PRODUCTS

TransMAX



TransMAX is a peptide-based solution to enhance uptake of adenovirus and increase transgene expression. New TransMAX enhances adenovirus transduction by specifically bridging the adenovirus surface to the cell membrane. Previously hard-to-transduce and even refractory cell types are now accessible to adenovirus transduction.

GFP Transduction

without



with TransMAX



NIH-3T3 cells were transduced with an adenovirus Serotype 5 expressing EGFP with a MOI of 200 with and without pre-incubation of TransMAX.

TransMAX:

- 20 to 50-fold enhanced expression
- Access to nearly all cell types
- Reduced adenovirus amounts
- No cell toxicity

Order Form

Complete this form and return it to: **SIRION BIOTECH GmbH, Am Klopferspitz 19, 82152 Martinsried.**
 Alternatively send the form by Fax: **0049 (0) 89 700 961 998**

Billing address

Delivery address

(only if different from billing address)

Name/Company

Name/Company

Street

Street

ZIP Code, City

ZIP Code, City

State

State

Tel/Fax

Tel/Fax

Email (for receiving an order confirmation)

Order

Order No.	Description	Amount	Price	Total
750-10	TransMAX Adenovirus Transduction Enhancer Solution, 50 Transductions		295 EUR	
750-20	TransMAX Adenovirus Transduction Enhancer Solution, 150 Transductions		695 EUR	
			Sub total:	
			Customer ID if known:	
			Freight charge:	
			TOTAL:	

Herewith I order the product(s) from the list above and agree with the general terms of SIRION BIOTECH GmbH (<http://www.sirion-biotech.de/COM/legal/agb.htm>).

Signature

Freight charge

Europe	29.- EUR
Overseas	59.- EUR

Please note: For orders from Germany the MWST of 19% will be added to the TOTAL.